



Standard glass type pH sensors



HI 1002/3

0-14 pH, double Teflon® junction, polymer electrolyte, -5 to 80°C, 6 bars (87 PSI) pressure, polypropylene body, BNC connector & 3 m cable.



HI 1002/5

0-14 pH, double Teflon® junction, polymer electrolyte, -5 to 80°C, 6 bars (87 PSI) pressure, polypropylene body, BNC connector & 5 m cable.



HI 1003/3

0-14 pH, double Teflon® junction, polymer electrolyte, matching pin, -5 to 80°C, 6 bars (87 PSI) pressure, polypropylene body, BNC connector & 3 m cable.



HI 1003/5

0-14 pH, double Teflon® junction, polymer electrolyte, matching pin, -5 to 80°C, 6 bars (87 PSI) pressure, polypropylene body, BNC connector & 5 m cable.



HI 1001

0-14 pH, double Teflon® junction, polymer electrolyte, -5 to 80°C, 6 bars (87 PSI) pressure, polypropylene body, BNC connector & 3 m cable.



HI 101

0-14 pH, double Teflon® junction, polymer electrolyte, -5 to 100°C, 6 bars (87 PSI) pressure, PVDF body & BNC connector.



HI 2910B/5

Double cloth junction, gel electrolyte, -5 to 80°C, 3 bars (43.5 PSI) pressure, Ultem® body, BNC connector & 5 m cable.

Platinum type ORP sensors



HI 2002/3

± 2000 mV, double Teflon® junction, polymer electrolyte, -5 to 80°C, 6 bars (87 PSI) pressure, polypropylene body, BNC connector & 3 m cable.



HI 2002/5

± 2000 mV, double Teflon® junction, polymer electrolyte, -5 to 80°C, 6 bars (87 PSI) pressure, polypropylene body, BNC connector & 5 m cable.



HI 2003/3

± 2000 mV, double Teflon® junction, polymer electrolyte, matching pin, -5 to 80°C, 6 bars (87 PSI) pressure, polypropylene body, BNC connector & 3 m cable.



HI 2003/5

± 2000 mV, double Teflon® junction, polymer electrolyte, matching pin, -5 to 80°C, 6 bars (87 PSI) pressure, polypropylene body, BNC connector & 5 m cable.



HI 2001

± 2000 mV, double Teflon® junction, polymer electrolyte, -5 to 80°C, 6 bars (87 PSI) pressure, polypropylene body, BNC connector & 3 m cable.



HI 2930B/5

Double cloth junction, gel electrolyte, -5 to 80°C, 3 bars (43.5 PSI) pressure, Ultem® body, BNC connector & 5 m cable.