

pH 1

○ pH

MAIN FEATURES

IP67 watertight portable pH meter. Resolution 0.1 pH with microprocessor and fixed sensor. Push-button calibration on 1 or 2 points (choosing between pH 4, 7, 10). Wide LCD display with indication of the buffers used for calibration and measurement stability. **Automatic temperature compensation.** Auto power off function. **Self-diagnosis messages: low batteries, calibration errors.** Power supply with 2 x 1.5V AAA batteries included, compatible with rechargeable batteries too. 6-month warranty on electronic parts.

Measuring range: pH 0,0...14,0 (0,1 pH). Working temperature: 0..50°C

The KIT includes: standard buffer solutions sachets pH 7 and pH 4 , 25 ml, sachets storage solution, calibration cap, lanyard, tissue, quick manual and cardboard case.



Name	Description	p/n
pH 1 KIT	pH 1 Kit complete	50014013
pH 1 ECOPACK	pH 1, batteries 1,5 V AAA, lanyard, tissue	50014003

TESTER

COND 1

○ Conductivity - TDS

MAIN FEATURES

IP67 watertight portable conductivity meter for Conductivity and TDS measurement. Indication of the buffers used for calibration and measurement stability. Auto power off function. Self-diagnosis messages. Power supply with 2 x 1.5V AAA batteries included, compatible with rechargeable batteries too. 6-month warranty on electronic parts.

Measuring range: COND 0,01 µS...200 mS; TDS 0,01 ppm...199,9 ppt. Working temperature: 0..50°C

The KIT includes: standard buffer solutions sachets 1413 µS e 12,88 mS 25 ml, lanyard, tissue, quick manual and cardboard case.



Name	Description	p/n
COND 1 KIT	COND 1 Kit complete	50014163
COND 1 ECOPACK	COND 1, batteries 1,5 V AAA, lanyard, quick manual	50014153